A Scalable Artificial Intelligence Framework for Rapid EGFR Mutation Screening in Lung Cancer

Neeraj Kumar^{1,*} Swaraj Nanda¹ Siddharth Singi¹ Gabriele Campanella² Thomas Fuchs² Chad Vanderbilt^{1,*}


¹Memorial Sloan Kettering Cancer Center, New York, NY 10065

²Icahn School of Medicine at Mount Sinai, New York, NY 10029

*Corresponding authors: kumarn6@mskcc.org (N.K.), vanderbc@mskcc.org (C.V.)







Memorial Sloan Kettering Cancer Center_{TM}

We present two Al approaches for EGFR mutation screening in lung

achieve 0.90 AUC on independent testing with 2-minute inference


time per slide. The in-real-time pipeline processes slides within 20-

cancer using H&E-stained whole slide images. Our methods

45 minutes of scanning, compared to conventional molecular



Whole Slide Image




Acknowledgements

We thank the Warren Alpert foundation and the department leadership for their support in this work.

Contact

- Chad Vanderbilt (vanderbc@mskcc.org)
- Neeraj Kumar (kumarn6@mskcc.org)



Approach 1

Two-step process: Pre-trained vision transformer (frozen) + MIL

Note: Bold AUC values indicate better performance


Approach 2

End-to-end training: Joint optimization of encoder and MIL

Dataset	Approach 1			Approach 2		
	AUC	Sens	Spec	AUC	Sens	Spec
Model Selection (n=260)	0.96	0.83	0.92	0.90	0.88	0.93
Independent Testing (n=6300)	0.89	0.72	0.96	0.90	0.84	0.95
TCGA-LUAD (n=519)	0.78	0.60	0.81	0.86	0.78	0.74
IRT (n=1000)	0.80	0.65	0.80	0.83	0.62	0.83

testing that takes 2-18 days.

EGFR Testing Methods: Turnaround Time



Note:

The AI-based method processes all slides within 30-45 minutes (median: 20 minutes), compared to days for conventional methods.